

SUB1-SUB1 homolog (S. cerevisiae) Gene

Reference: PTXBC010537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SUB1 homolog, transcriptional regulator; SUB1 homolog; P15; PC4; p14; activated RNA polymerase II transcriptional coactivator p15; activated RNA polymerase II transcription cofactor 4; positive cofactor 4		

For research use only.

GENE DESCRIPTION	SUB1 homolog (S. cerevisiae)
GENE ID	10923
ORIGIN SPECIES	Human
ACCESSIONS	BC010537
VECTOR	pDONR223
NCBI SYMBOL	SUB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctaaatcaaaaggaactgtttcttcaagctcttctggcagtgattctgacagtgaggtgacaaaaagttaaaggagaaaaagcaagttgtccagaaaaacctgtaagaacaaaaagacaggtgagacttcgagagccctgtcatcttctaaacagacagcagcagcagagatgat
 aacatgtttcagattgggaaatgaggtacgttagtgttcgcgattttaaaggcaaatgctaatgatatattagagaataatggatggatcctgaagtgaaatgaaccaggaagaaaagggtatttcttaaatccagaacaatggagccagctgaaggacagattctgacattgatgatgca
 gtaagaaaactgtaa