

## ZNF167-zinc finger protein 167 Gene

Reference: PTXBC010555

|                |                                                                                                                                                                                                           |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ZNF167; ZFP; ZNF448; ZNF64; ZSCAN39; zinc finger protein with KRAB and SCAN domains 7; zinc finger protein 167; zinc finger protein 448; zinc finger protein 64; zinc finger with KRAB and SCAN domains 7 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger protein 167                                                                           |
| <b>GENE ID</b>          | 55888                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC010555                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ZNF167                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaccactgcaggcagggaatttaggcctcatcccaggagcactgcttccagaagcaagagggcgccctgactgtgaagcaggagccagcaaacagacctgggggcagggcagcagctctccagaagaactatcctctgtctgcgaatctccggctacactcaggcaatt  
gtgttaccacagatgtctggccgcagggaagcattgagccgcttccggagctctgccgtggtgctcatgccagaggtgcacaccaaggagcagatcctggagctgctgtgctgagcagttcctgagcatctccctgggagctccgacctgggtgcagctgcatccctg  
agagtgtgagggagctgtggtgtggtgaggttccagagacacctcagtgatcagaggaggtttcagccccagcacagaacaggaatgcatgttggagagacaacagctctgggtacaacaagggaatctcctctacacctacccctcagtggggctcagccctggagcc  
cacctggagcctcttatgacccagggacacaccacctccagtggggagctcgtcaatgtacttctcagttctcctcctcaagtggggaactcaggagaccaagcaggggcaactgtacttggatgtcaggccccaggatactgtggcatatgaggacctatctgtagacta  
cactcagaagaatggaaaagtcacactcagtcagagagccctgcagtgaacatgatgccagaaatcaccatagcatggcctccttgggctggagtacaatggcgtga