

FAM107A-family with sequence similarity 107, member A Gene

Reference: PTXBC010561

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM107A; DRR1; TU3A; down-regulated in renal cell carcinoma 1; family with sequence similarity 107 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 107, member A
GENE ID	11170
ORIGIN SPECIES	Human
ACCESSIONS	BC010561
VECTOR	pDONR223
NCBI SYMBOL	FAM107A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtactcggagatccagagggagcgggcagacattggggccctgatggccggccagaatacagagagtggaatcgggagctcatcaagcccaagaagctgctgaacccctgaaggcctctcggagtcaccaggagctccaccgggagctgctcatgaaccacagaaggggccc
ttggtgtggacagcaagccagagctgcagcgtgtcctagagacccggcggcgaaccagctcatcaagaagaagaaggaggagctggaagccaagcggctgcagtgccctttgagcaggagctgctgagacggcagcagagggctgaaccagctggaaaaaccacagagaag
gaagaggatcacgccccgagtttataaagtcagggaacacctgcggagaattgccacactgaccagcgaagagagagagctgtag