



Reference: PTXBC010569

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16DUPq36.3; DUP7q36.3; PACAP-R-3; PACAP-R3; VIP-R-2; VPAC2; VPAC2R; VPCAP2R; vasoactive intestinal polypeptide receptor 2; PACAP type III receptor; VIP and PACAP receptor 2; helodermin-preferring VIP receptor; pituitary adenylate cyclase-activating polyp		

For research use only.

GENE DESCRIPTION	vasoactive intestinal peptide receptor 2
GENE ID	7434
ORIGIN SPECIES	Human
ACCESSIONS	BC010569
VECTOR	pDONR223
NCBI SYMBOL	VIPR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]