

IGF1R-insulin-like growth factor 1 receptor Gene

Reference: PTXBC010607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	soluble IGF1R variant 2; soluble IGF1R variant 1; CD221; IGFIR; IGFR; JTK13; insulin-like growth factor 1 receptor; IGF-I receptor; insulin like growth factor 1 receptor		

For research use only.

GENE DESCRIPTION	insulin-like growth factor 1 receptor
GENE ID	3480
ORIGIN SPECIES	Human
ACCESSIONS	BC010607
VECTOR	pDONR223
NCBI SYMBOL	IGF1R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaacagccgaggtgttgagcccagcagtgcatggcaccgttcggcatctgcttgattggtctggctgccgtcattgtcagcagtgccatggacatgggaagacttgactgcacagccaatggtttcatgatgattacagcatacacagtgatcacataa