

C14orf126-chromosome 14 open reading frame 126 Gene

Reference: PTXBC010618

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf126; D-tyrosyl-tRNA deacylase 2 (putative)		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 126
GENE ID	112487
ORIGIN SPECIES	Human
ACCESSIONS	BC010618
VECTOR	pDONR223
NCBI SYMBOL	C14orf126
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgagggtagccggattctcaggcccgccgctctacagcagtgccgcagcccgctgcaaatcgccagccgagtgaggacgtcgccgcccagtgagggtgaggtccaaagaggactggtgatctacgtgtgtttcaaggaggctgataaagaacttcttccaaaat
gggttaatacactgttaaatgtgaaattaagtgagacagaaaatggcaagcatgtctctatattggatctacctggcaacattcttattatccctcaagctaccctggagggaagactaaaaggagaacatgcaatatcactcttaactctgaaaaagaagggtttgaactttacttcaattgt
gactctatgtgaaaaagaagtagctgctaatagcaagtgctgaagctagggtttagtggaacatggcacttatgggaacaggcaggtgttaagctggacaccaacggaccattcacacacttaattgagtttga