



Reference: PTXBC010645

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGS23; NANTA3; SAT3; SIAT4; SIAT4C; ST3GalIV; STZ; CMP-N-acetylneuraminase-beta-galactosamide-alpha-2,3-sialyltransferase 4; SAT-3; SIAT4-C; ST-4; ST3GalA.2; alpha 2,3-ST 4; alpha 2,3-sialyltransferase IV; alpha-3-N-acetylneuraminyltransferase; beta-galac		

For research use only.

GENE DESCRIPTION	ST3 beta-galactoside alpha-2,3-sialyltransferase 4
GENE ID	6484
ORIGIN SPECIES	Human
ACCESSIONS	BC010645
VECTOR	pDONR223
NCBI SYMBOL	ST3GAL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]