

SNAP25-synaptosomal-associated protein, 25kDa Gene

Reference: PTXBC010647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMS18; RIC-4; RIC4; SEC9; SNAP-25; SUP; bA416N4.2; dJ1068F16.2; synaptosomal-associated protein 25; resistance to inhibitors of cholinesterase 4 homolog; synaptosomal-associated protein, 25kDa; synaptosome associated protein 25kDa; synaptosome associated		

For research use only.

GENE DESCRIPTION	synaptosomal-associated protein, 25kDa
GENE ID	6616
ORIGIN SPECIES	Human
ACCESSIONS	BC010647
VECTOR	pDONR223
NCBI SYMBOL	SNAP25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaagacgcagacatgcgcaatgagctggaggagatgcagcgaagggtgaccagtggctgatgagtcgctggaagcaccgctgctgctgcaactggtgaagagagtaaatgctggtatcaggacttgggtatggtgatgaacaaggagaacaactcgatcgtgc
gaagaaggcatgaaccatcaaccaagacatgaaggagctgagaaaaatttaaaagatttagggaaatgctgtggcctttcatatgtcctgttaacaagcttaaatcaatgagctttacaaaaagcctggggcaataatcagacggagtggtggccagccgctgctgtagtg
gacgaacgggagcagatggccatcagtggtgcttcccgagggttaacaaatgatgccgagaaaaatgaatggatgaaacctagagcaggtgagcggcatcatcggaacctccgtcacatggccctggatgggcaatgagatcgatacacagaatgccagatgcacagg
atcatggagaaggctgattccaacaaaccagaattgatgagccaacacgtgcaacaaagatgctgggaagtgttaa