

LRRFIP1-leucine rich repeat (in FLII) interacting protein 1 Gene

Reference: PTXBC010662

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLAP-1; FLAP1; FLIIAP1; GCF-2; GCF2; HUF1-1; TRIP; leucine-rich repeat flightless-interacting protein 1; GC-binding factor 2; LRR FLII-interacting protein 1; NEDD8-conjugating enzyme; TAR RNA-interacting protein; leucine rich repeat (in FLII) interacting		

For research use only.

GENE DESCRIPTION	leucine rich repeat (in FLII) interacting protein 1
GENE ID	9208
ORIGIN SPECIES	Human
ACCESSIONS	BC010662
VECTOR	pDONR223
NCBI SYMBOL	LRRFIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatgtaatggacctacaaggatgccacagacagatcagcgacctcaaatftaaactgcataatctgagcaagagataactgcattagaacaaatgtaataaggttagagagtcaagatcacgttacaaatcagcggtgaaaatgcagaaaaatagaagatgaacttaagc
agaaaaacggaaactcacaagagagctccgctctgattggataaaacagaagagctcgaggtgagcaacggccacttagtgaagcgtctgaaaaatgaaagcaaatcggagtgcactctgtccagcagtaa