

NDUFB11-NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 11, 17.3kDa Gene

Reference: PTXBC010665

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-ESSS; ESSS; NP17.3; Np15; P17.3; NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 11, mitochondrial; NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 11, 17.3kDa; NADH-ubiquinone oxidoreductase ESSS subunit; complex I NP17.3 subunit; complex		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 11, 17.3kDa
GENE ID	54539
ORIGIN SPECIES	Human
ACCESSIONS	BC010665
VECTOR	pDONR223
NCBI SYMBOL	NDUFB11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctgggctgttgggttgagcgtcgccgtcttttggcggcagcggcgacgcgaggggtccggccgcccgcgtccgctgggaatctagcttctccaggaactgtgttcgccccgtccgctgtggcgggaaagcggccccagaccgaccacaccgtggcaagaggaccca
gaaccggaggacgaaactgtatgagaagaaccagactcccatggttatgacaaggacccgttttggacgtctgggaacatgcgactgtcttctcttggcgtctccatcatctgtgtcttggcagcaccttttggcctaictgcctgactacaggatgaagagtgggtcccgccgga
agctgagaggcttgaataaccgagagccaatggcctcccatcatggaatccaactgttcgacccagcaagatccagctgccagaggatgagtga