

MLX-MAX-like protein X Gene

Reference: PTXBC010689

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MLX, MAX dimerization protein; MAD7; MXD7; TCFL4; TF4; bHLHd13; max-like protein X; BigMax protein; MAX-like bHLHZIP protein; class D basic helix-loop-helix protein 13; transcription factor-like protein 4		

For research use only.

GENE DESCRIPTION	MAX-like protein X
GENE ID	6945
ORIGIN SPECIES	Human
ACCESSIONS	BC010689
VECTOR	pDONR223
NCBI SYMBOL	MLX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggagccggcgccctcccagggacccttgggtcaagggtgagtagtgcctacagcgacaacagcctggacccgggctttttagaagacccgcaaggggagtgtagtgcagagcgaatagcatcggtccaccagtgccctctctgtcccccaacacagatgatgagga
cagtgattaccaccaggaggcctacaaggagtcctacaagaccggcgccggcgccgcacacacacacagcctgagcagaagaggaggagccatcaagagaggtatgatgacctcagaccatcgccccacttgccagcagcaggacttctccattggtcccaaaagctcagcaa
agccatcgttctacaaaagaccattgactacattcagttttgcacaaggagaagaaaaagcaggaggagggtgtccacgttacgcaaggatgtcaccgccctaaagatcatgaaagtgaactatgagcagattgtgaaggcacaccaggacaacccccatgaaggggaggaccag
gtctctgaccagggtcaagttcaacgtgtttcaaggcatcatggattccctgttccagtccttcaatgccctccatcagtgccagcttcagagagctgtcagcatgtgtcttcagctggatcagggagcactgtaagcctcagaccctgcgggagattgtattggcgtcctgcaccaattgaaa
aaccagctttactga