

LOC404266-hypothetical LOC404266 Gene

Reference: PTXBC010732

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOXB cluster antisense RNA 3		

For research use only.

GENE DESCRIPTION	hypothetical LOC404266
GENE ID	404266
ORIGIN SPECIES	Human
ACCESSIONS	BC010732
VECTOR	pDONR223
NCBI SYMBOL	LOC404266
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagtctgccgggaaccagcgatccgaccagcggagaaggttcaggctgccggcgccggcgagagagcgggaagagagctcggaggagccccggcgctggcgtggtcaggctccgagagcggccgggatcgccacaccggcctgtaa