



Reference: PTXBC010863

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTPMEG2; tyrosine-protein phosphatase non-receptor type 9; PTPase-MEG2; protein-tyrosine phosphatase MEG2; protein tyrosine phosphatase, non-receptor type 9		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 9
GENE ID	5780
ORIGIN SPECIES	Human
ACCESSIONS	BC010863
VECTOR	pDONR223
NCBI SYMBOL	PTPN9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr