

EIF4EBP3-eukaryotic translation initiation factor 4E binding protein 3 Gene

Reference: PTXBC010881

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4E-BP3; 4EBP3; eukaryotic translation initiation factor 4E-binding protein 3; eIF4E-binding protein 3; eukaryotic initiation factor 4E-binding protein 3; eukaryotic translation initiation factor 4E binding protein 3		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 4E binding protein 3
GENE ID	8637
ORIGIN SPECIES	Human
ACCESSIONS	BC010881
VECTOR	pDONR223
NCBI SYMBOL	EIF4EBP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaacgtccactagctgccgattccggggggccgggaccagctgccgactgtctacagaccacgccggggggcagctatacgcactacccccggaggcaccaggatcatctacgacgaaagtctctgtggagtgcagaactcaccattgccggacacccccctgc
tgcctccctcagattccgggggtcacaactctccaacagccccctctccaagctggaggagctgaaggagcaggagacagaggagagatacccgatgacgcacaattgaaatggacatctaa