

CCDC53-coiled-coil domain containing 53 Gene

Reference: PTXBC010889

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WASH complex subunit CCDC53; CCDC53; CGI-116; coiled-coil domain containing 53; coiled-coil domain-containing protein 53; WASH complex subunit 3		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 53
GENE ID	51019
ORIGIN SPECIES	Human
ACCESSIONS	BC010889
VECTOR	pDONR223
NCBI SYMBOL	CCDC53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgaggacgggtctctcatgggtcaggcatagacctgaccaaggtgccagctattcaacagaaagaacgggtgctttttaaaccaattgtggtgcacactgtacagttcctcaaccgttttctacagttgtgaggagaaactggcagaccttctcgtatccaacaatt
gaaacaactctcaatatttttagatgcaaaagttgtcatctatccaggcctagatggtcacagttgaagtatctctttaaattgacacagtggtcacaaatggagcacatcctgaagccacttcagagcaaccacagcagagcagtagacacagactctggactacaggaaagtgaagtatcag
cagaaaaatcttcaactgtagccaaaggtccaagatagccagatattcaaaatggttcaagtggtgtaccagtgatggcaataagaacaaaatgatacagaaggactagaccagatcttctgagaggccagatgctccagtgctgatggcgaaagtgaagaaactgtagaagaa
agttcagatagcgaatctcttttagtgattaa