

COMMD9-COMM domain containing 9 Gene

Reference: PTXBC010892

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC166; COMM domain-containing protein 9; COMM domain containing 9		
For research use only.			

GENE DESCRIPTION	COMM domain containing 9
GENE ID	29099
ORIGIN SPECIES	Human
ACCESSIONS	BC010892
VECTOR	pDONR223
NCBI SYMBOL	COMMD9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgccctgacagcggagcattttgcagcactccagagcctgctcaaggcctctcgaagatgtttgcagacagctgtgcaagaaagctttccagttcagccctggcttgaaaaaaccttggatgttacatgttccagctgtctgtgacccaggaggagggcagaggaaactgctcc
aggctctgcaccgcctcactagctgggtgcatcctgtgacctgtcctctgccgaggaattctggctctcttccagaaaaattccacaaaacctcaaaaacctgctgacaaagatcatctagaacatgtgtctacttggagaaccgaagccaggcgaatcagatctctgccacgcct
ggctgatctggactggagatggatatacaaacctctcagacagcatcagccgcatggcctgccccacctgctgctccagatgaagatccaagaagatccagcctatgaggagacaaacctccatctcagctgtcaccgtggagctgagcaagaacactggacacatgttag
atggcctggccgcgacccagacaaactctctgcctggccagtaaatga