

## DEPDC1B-DEP domain containing 1B Gene

Reference: PTXBC010904

|                |                                                                                                                                                                                                      |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | BRCC3; XTP1; DEP domain-containing protein 1B; HBV X-transactivated gene 8 protein; HBV XAg-transactivated protein 8; HBxAg transactivated protein 1; breast cancer cell 3; DEP domain containing 1B |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | DEP domain containing 1B                                                                          |
| <b>GENE ID</b>          | 55789                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC010904                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | DEPDC1B                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgccacccctgtgtgatggcttggtacccgaacactgatgttcagacatttcccggttcacatctgtgtccaaggatgaagtggaacttgatgagttattagctctagattgtaacgtttctgatggacaattaccaggaaattctgaaagtcctttggccttcagacctctatagaggag  
cgtgtggtctatctacgaagagtcagataaaataccaggagctgatgtgatcactttatctgctccatcatttgcgtcaaatagtcagaggaatttgaatatcaagatcatatggctctcaggaaacctctggcagcctgttggagggaagtataacagatgccaaactctccaaca  
aagagaaaaagaagaactgaagcagtttcagaatctctcctgaagtcctatcaagaacgatttcctacaccagaagtgacgacattctgttctgaaaaacccaaccgaaccacagctgtaattgtggcactaaagaagcctttcaaccattcaagaactagaagtttctgaat  
gtaa