

MED9-mediator complex subunit 9 Gene

Reference: PTXBC010906

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MED25; mediator of RNA polymerase II transcription subunit 9; mediator of RNA polymerase II transcription, subunit 9 homolog; mediator subunit 25; mediator complex subunit 9		

For research use only.

GENE DESCRIPTION	mediator complex subunit 9
GENE ID	55090
ORIGIN SPECIES	Human
ACCESSIONS	BC010906
VECTOR	pDONR223
NCBI SYMBOL	MED9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctctgctgggtggcagccggcgacaggcgaggatgtattgccccaacgtccgaccagccgctgacctgacaccaagccgctgcccctcctcagccgcccgggtccctgcgcctcaaccgcagcagtcgcccggcgccacggcctcagtcacctgcccgcgaggagggaagagaactactcctttttacctttggttcacacatcatcaaatggtaa