

C1orf123-chromosome 1 open reading frame 123 Gene

Reference: PTXBC010908

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0587 protein C1orf123; chromosome 1 open reading frame 123		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 123
GENE ID	54987
ORIGIN SPECIES	Human
ACCESSIONS	BC010908
VECTOR	pDONR223
NCBI SYMBOL	C1orf123
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaaatcgcgctgcaactcaaagccagctggagaacatcaccaactccggccgtggcgaggactccggtgtacctgaagatgaaatgtggcaactgtggtgagattcggacaagtgccagctacatccgctgatggacagtgtggcactgaagggggccgtggc
agtgctccatggtccagaagtgcagctgtgtgcaagagaaaattccatcgagattttaagcagcaccatcaagccttacaatgctgaagacaatgagaactcaagacaatagtgagtttgagtgccggggcctgaaccagttgattccagccagggctgggttgcctgaaggt
gtggagtcagggacagccttcagtgcatttaattcgcaggagaaggactggactgactatgataaaaggccaggagctgtgggaatctatgaggtcaccaccagtttgaagtgctga