

COX17-COX17 cytochrome c oxidase assembly homolog (S. cerevisiae) Gene

Reference: PTXBC010933

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX17, cytochrome c oxidase copper chaperone; COX17 cytochrome c oxidase assembly homolog; cytochrome c oxidase copper chaperone; cytochrome c oxidase 17 copper chaperone; cytochrome c oxidase assembly homolog 17; human homolog of yeast mitochondrial copp		

For research use only.

GENE DESCRIPTION	COX17 cytochrome c oxidase assembly homolog (S. cerevisiae)
GENE ID	10063
ORIGIN SPECIES	Human
ACCESSIONS	BC010933
VECTOR	pDONR223
NCBI SYMBOL	COX17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgggtctgtgactcaaacctgccccgcctgagtcagggagaagaagccgctgaagccctgctgcgttgcctggagaccaagaaggcgcgcgatgcgtgtatcatcagagaaggagaagaacactgtggacatctaattgagggccacaaggaatgcatgagagcccta
ggatttaaaatatga