

PI3-peptidase inhibitor 3, skin-derived Gene

Reference: PTXBC010952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ESI; SKALP; WAP3; WFDC14; cementoin; PI-3; WAP four-disulfide core domain 14; WAP four-disulfide core domain protein 14; elastase-specific inhibitor; peptidase inhibitor 3, skin-derived; pre-elafin; protease inhibitor 3, skin-derived (SKALP); protease inh		

For research use only.

GENE DESCRIPTION	peptidase inhibitor 3, skin-derived
GENE ID	5266
ORIGIN SPECIES	Human
ACCESSIONS	BC010952
VECTOR	pDONR223
NCBI SYMBOL	PI3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggccagcagcttcttgatcgtggtggtgttcctcctcgtgggacgtggtctagaggcagctgtcacgggagttctgttaaagggtcaagacactgtcaaggccgtgttcattcaatggacaagatcccgtaaaggacaagttcagttaaagggtcaagataaagtcaggcgc
aagagccagtcagggtccagctccactaagcctggtcctgcccattatcttgatccgtgctgcatgttgaatccccctaaccgctgcttgaagatactgactgccagggaatcaagaagtgtgtgaaggctcttgcggatggcctgtttcgttccccagtga