

PFDN4-prefoldin subunit 4 Gene

Reference: PTXBC010953

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PFD4; prefoldin subunit 4; prefoldin 4; protein C-1		
For research use only.			

GENE DESCRIPTION	prefoldin subunit 4
GENE ID	5203
ORIGIN SPECIES	Human
ACCESSIONS	BC010953
VECTOR	pDONR223
NCBI SYMBOL	PFDN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggccacatgaagaaggcggctgcagaagatgtcaatgttactttcgaagatcaacaaaagataaacaatttcacgggaatacaagtagaatcacagagctgaagggaagaatagaagtaaaaaagaacaaactccaaacctagaagatgcttgtgatgacatcatgcttgcagatgatgtcttaataatgatacccttatcaaaattggtgatgtcttcattagccattctcaagaagaacgcaagaatgttagaagaagcaagaaaaattgcaagaagaattgacgcttagaatccagagtggaatcaattcagcagtgtagcagattgaaagtcagttgtatgcaaaatcgggagcaacataaaccttgaagctgatgaaagttaa