

ZFAND5-zinc finger, AN1-type domain 5 Gene

Reference: PTXBC011018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZA20D2; ZFAND5A; ZNF216; AN1-type zinc finger protein 5; zinc finger A20 domain-containing protein 2; zinc finger protein 216; zinc finger, A20 domain containing 2; zinc finger, AN1-type domain 5; zinc finger AN1-type containing 5		

For research use only.

GENE DESCRIPTION	zinc finger, AN1-type domain 5
GENE ID	7763
ORIGIN SPECIES	Human
ACCESSIONS	BC011018
VECTOR	pDONR223
NCBI SYMBOL	ZFAND5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcaggagactaaccagacccggggcccatgctgtgtgacacaggtgtggctttatggaaatcctaggacaaatggaatgtgtcagttgtctacaaagaacatctcagaggcagcaaaatagtgccagaatgagcccaatggggacagctagtgttccaacagctctacctc
agattctgcatctgtacagagagcagacactagcttaacaactgtgaaggtgtctgtgacacacatctgaaaaatcaagaatgtgcctgtgctgcctgttaactcagcaaatgacagaatgagcatttcaagagaggacaaaataactacccgaaaacagaggtgtcagag
ccagttgtcactcagccagtcctcagcttctcagccagctctcagagtgaaagaaagctctgaattgcccaaccaagaaaaacagatgttcatgtgcagaagaaagttgtgttactcaggggttgactgccgatgtggaattgtttgtggacttcaccgttactctgacaagc
acaactgtccgtatgattcaaaagcagaagctgcagcaaaatcagaaaagagaatccagttgtgtgctgaaaaaattcagagaataataa