

POFUT2-protein O-fucosyltransferase 2 Gene

Reference: PTXBC011044

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FUT13; GDP-fucose protein O-fucosyltransferase 2; O-FucT-2; peptide-O-fucosyltransferase 2; protein O-fucosyltransferase 2		

For research use only.

GENE DESCRIPTION	protein O-fucosyltransferase 2
GENE ID	23275
ORIGIN SPECIES	Human
ACCESSIONS	BC011044
VECTOR	pDONR223
NCBI SYMBOL	POFUT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgttgcaggacactgcgggaggtgggagacgagttcaggagcagacatcctaactccaggacgacgacagagatcccttcaggaggactggatgaagatgaaggtcaagctgggtccgcgtagggggccctacctgggagtcacctgagaagaaaagattc
atctgggtcacagacaggatgtaccagctctggaaggccgtgaggaaagatccgagcctcatgaagaccaccgctggacaaggtgttggccacagatgccgtcagaaggaataatgaagagctaaaaaagctgttaccgagatggtgaggtttgaaccacgtgggag
gagctggagctctacaaggacggagcggttgcgattatgaccagtggtatctgcgcacagccagttcatga