

TYROBP-TYRO protein tyrosine kinase binding protein Gene

Reference: PTXBC011175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KARAP; PLOSL; TYRO protein tyrosine kinase-binding protein; DNAX-activation protein 12; KAR-associated protein; killer-activating receptor-associated protein; TYRO protein tyrosine kinase binding protein		

For research use only.

GENE DESCRIPTION	TYRO protein tyrosine kinase binding protein
GENE ID	7305
ORIGIN SPECIES	Human
ACCESSIONS	BC011175
VECTOR	pDONR223
NCBI SYMBOL	TYROBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggggaactgaacctgcagcaggtctctgctctcctcctgctggtgtaagtgtctcctcctgtccaggccagggccagagcattgagttgctctacggtgagccggggtgctggcagggatcgtgatgggagacctggtgctgacagtgtcattgccctggccgtgtacttctgggcccgtggtccctcggggcgaggggctgcggaggcagcgaccggaaacagcgtatcactgagaccgagtcgccttatcaggagctccagggtcagagtcggatgtctacagcgacctcaacacagagggcgtattacaatga