



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PPP2R3B-protein phosphatase 2 (formerly 2A), regulatory subunit B'', beta Gene

Reference: PTXBC011180

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

NYREN8; PPP2R3L; PPP2R3LY; PR48; serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta; NY-REN-8 antigen; PP2A, subunit B, PR48 isoform; protein phosphatase 2 (formerly 2A), regulatory subunit B", beta; protein phosphatase 2, regula

For research use only.

GENE DESCRIPTION	protein phosphatase 2 (formerly 2A), regulatory subunit B", beta
GENE ID	28227
ORIGIN SPECIES	Human
ACCESSIONS	BC011180
VECTOR	pDONR223
NCBI SYMBOL	PPP2R3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatagacagagatcttcaggagcagtcacacaggcagaaaagtgcgaaggagggaagatcagctatgccactttgtctgttttgcctctgcgaagacaataaacaccgaccagcatcagctactgttctgcctgcatacctgacggggacgcgcctcttccatg
ttcgagctcgagctactctacgaggagcagtgccgaagctcagcactatgcccactgaggccctgcccttcaggaactccttgcagatgctgacctgtgtcaagccgaggactgaagggaatcacctctgcgaagcctgaagcgtctcaagctggccaacgtcttcttcgaca
ctctcttcaactcagagaagctacctgcaccagcagagaagcagcatctccctgtcaggcagcgtgtgcagcggcgcccccagctcctggactgcgaagaaagtacgcggccgagagatcagcctctgtgtgcggaggagagcctggcgagagcccttggagagagcgtgttcg
agcccgccagccactcctgtgcagagaagctgtagtgcctgcctccctgcggccagagagcccttctgcagcgcctcactcagctggcgccgtgtgaacctgtacagctacatgcacatgcggagacgagcactgcagcctgcagcctctga

More info: Please contact us for any question: contact@proteogenix.fr