



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HIST1H2BN-histone cluster 1, H2bn Gene

Reference: PTXBC011372

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2B/d; H2BFD; histone H2B type 1-N; H2B histone family, member D; histone 1, H2bn; histone H2B.d; histone cluster 1, H2bn; histone cluster 1 H2B family member n		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2bn
GENE ID	8341
ORIGIN SPECIES	Human
ACCESSIONS	BC011372
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2BN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccgagccctcaaaagctcgctctgccccgaagaagagctccaagaaggcagtgacaaaagcccgagaagaagcagcgcaagaagcgcaagcgagccgcaaggagagctactccgtgtactgtgtacaaaggtgtgaagcaggtccaccgccgacccggtatctgtctcaag
ggcatgggcatgatgaactctctgtcaatgacatcttctgcgcgtatccgcggcgagctgctcccgctgagcggtacattacaacagcgcgtcgatcattacacagcgcgtcgctgtctgcaggaagctgctgcgaagcagcggtgctgcgaagcagcggtgctgcgagggcacaag
cgcttcacagctacacacagctccgaagaagaagaaaaagatctcacagaagaacattagaagaagccgatgatgattcagaagaatccaaagtgtctgaattaaagctgcagcaacattccagatatacatattatgtccagagattga