



Reference: PTXBC011379

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C15orf39; chromosome 15 open reading frame 39		
For research use only.			

GENE DESCRIPTION	chromosome 15 open reading frame 39
GENE ID	56905
ORIGIN SPECIES	Human
ACCESSIONS	BC011379
VECTOR	pDONR223
NCBI SYMBOL	C15orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr