

C7orf29-chromosome 7 open reading frame 29 Gene

Reference: PTXBC011406

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf29; ZBED6 C-terminal-like protein; ZBED6 C-terminal like		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 29
GENE ID	113763
ORIGIN SPECIES	Human
ACCESSIONS	BC011406
VECTOR	pDONR223
NCBI SYMBOL	C7orf29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggtccgcgagggcagtgccgacagggcctctctgagccaggtgctgccccagctgcgctacctgcacatctctggagcaggttcacacactttcaggagcagagtgctcggggagagggcgagccatccagctggctgaggggttggccggcagctctgcaccgac
tgtcagctcaacaagctctctaccgcgagggagttgtgctggccacctgtgctgacctgtgtcaaggggaagattgagccatcctgcatggggccgaccgacatgaccactggaagcaagtcctctgttacaaggtgaaggagattagagtgctgaatactttgaactccca
agtccctgcaaagccccaggggtctgtgctggacccaccagggtagccaagctcgggtggaggggagaagccagggggagcctctgcagagcagcagccactctgggcttctgctggccagaggagaaaggcttctgagagcatggggctgtgctgctct
gagagaagcggggctgctgtccaccaagagccactggccagcatctgtcaagaagtatctgtgggagaatgagaccgttgagcccgatgacccctggttactgggagaagaagcgagaagcctggccacatctctgttaccctccacagagccttctga