

STMN4-stathmin-like 4 Gene

Reference: PTXBC011520

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RB3; stathmin-4; stathmin-like 4; stathmin-like protein B3; stathmin-like-protein RB3; stathmin 4		

For research use only.

GENE DESCRIPTION	stathmin-like 4
GENE ID	81551
ORIGIN SPECIES	Human
ACCESSIONS	BC011520
VECTOR	pDONR223
NCBI SYMBOL	STMN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgtctgcctacaaagagaagatgaaggagctccgctggtgtctctgtctctgctccgcatccctgaataagctgctacaaatgaaggctggtgtggagacagttaggaggaaggatgaaagccagcggaaagacagtgctgactggagagaaaga
agagctcaggcagacacgggtgacctgaattggtgctcatctccgacatggaagtcacgagctgaacaaatgcacctcgggccaatccttgaagtcacctgaagccaccctccttgatgggttcccgagttcaacgcctccctgccaagcggcgagaccatccctggaagaga
tccagaagaactagaagcggctgaggagcgaaggagtagcaggaagcggagctctgaaacacctagcagagaacgggaacatgagagagaggtgatccaaaggccattgaggaacaacaactcatcaagtggttaaggaaaaactggccagaagatggaatcca
acaaggagaacaggagccaccctcgcgcctatgttgaacggctgcaagagaaggacaagcagccgaggaggtgcgaaaaacaaggagctgaaggaaaggcctcaggttaa