



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD164-CD164 molecule, sialomucin Gene

Reference: PTXBC011522

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD164 molecule; CD164 antigen, sialomucin; DFNA66; MGC-24; MGC-24v; MUC-24; endolyn; sialomucin core protein 24; multi-glycosylated core protein 24		

For research use only.

GENE DESCRIPTION	CD164 molecule, sialomucin
GENE ID	8763
ORIGIN SPECIES	Human
ACCESSIONS	BC011522
VECTOR	pDONR223
NCBI SYMBOL	CD164
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcggtctccgcgtactcgtgtttggggcgcgcactcgtctggcggtcgtctgtgtgtgtctgcgcggagacaagaacgaccgaccggaacgtgacgactttagcgcgccattccaacgtaacctggcgccggtgacgtccctcccgctgtgaccactcgggcaccagaac
ctgtgaaggctgaaacagctcgtgttctctgtttaattgttagcgtgtgtaatactactcgttttgtagataaagttaaagatgagagctattgttcataactcaacagttagtgtgattgcaagtggggaaacgacagactctgtctcgtttccacggccactccagtgccaacagccaattctac
agctaaacccaacagttcagccctccctcttacaactccaagacagttactactcaggtacaaacaataaacactgtgactccaactcaacaactctgcgcaaggttactcctttgatgcagacttttcattggaggaggtctcctgtctcgtgtgcagcgtgtaatttctcttataaattct
gcaaatctaaagacacgaattaccactcgtctaa