

C14orf79-chromosome 14 open reading frame 79 Gene

Reference: PTXBC011537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C14orf79; chromosome 14 open reading frame 79		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 79
GENE ID	122616
ORIGIN SPECIES	Human
ACCESSIONS	BC011537
VECTOR	pDONR223
NCBI SYMBOL	C14orf79
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagattgtgacctcaagagcctgaaggactcctcactgtcagcagcttctgtctccagcattgcaaagccctgatccagaccaagctctcggggccgctggcagcaaacagggaggctgatgacatgcagccgcttctgaagacccctcatgcggaggtggccagacat
cactattccaaggaaaaggatgttctcctccagcaagctcaattgacactctttaatagcgacgtttgctaa