

ATP5S-ATP synthase, H⁺ transporting, mitochondrial F0 complex, subunit s (factor B)
Gene

Reference: PTXBC011549

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPW; HSU79253; ATP synthase subunit s, mitochondrial; ATP synthase coupling factor B, mitochondrial; ATP synthase coupling factor B-like 1; ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit s (factor B); ATP synthase-coupling factor B; mit		

For research use only.

GENE DESCRIPTION	ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit s (factor B)
GENE ID	27109
ORIGIN SPECIES	Human
ACCESSIONS	BC011549
VECTOR	pDONR223
NCBI SYMBOL	ATP5S
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgtgctgctgtgagcagcgactcactgtgcagatcaaatgatgtgtgtggaaaaattccagcagttgtgtggcgtaaagaaactccatggtcatgtgactccagatactctggggctggttgatgcagttgtaataagggtgattatgatcgcatcagggtgtggccctg
ataggcgccatccgagtggtgtgtgctgtgtggccatggtgcgtaccatggccaggagaggtggcagaaggactacaacccttcaacaggccctctggacaatacaagattcaggcgatcgacgccaccgactctgtatcatgagcattgattgatcacatggaacct
caaatattgtgttaa