

TMEM222-transmembrane protein 222 Gene

Reference: PTXBC011579

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf160; transmembrane protein 222		
For research use only.			

GENE DESCRIPTION	transmembrane protein 222
GENE ID	84065
ORIGIN SPECIES	Human
ACCESSIONS	BC011579
VECTOR	pDONR223
NCBI SYMBOL	TMEM222
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaagtggaggccgacggcgccgagacggacatgaagcaatatcaaggctccggcgctcgccatggatgtggaacggagtcgtccctactgcgtggtgtggagcccatcccggtgctcagtggtttttcccatcagcgccacatgggcatctgcacatcca
caggagtcattcgggacttcgcgccctactttgtctcagaggacaacatggcctttggaagcctgccaaagtactggaagttggacccctgctcaggtctatctagcgggcccacgcatgggacacggctgtgcacgacgacctgaggagtacaagcaccgcatgcacatctct
gctgtgacaactgccactcgcacgtggcattggccctgaatctgatgcgctacaacaacagcaccaactggaatatgtgacgctctgcttctctgctctacgggaagtacgtcagcggttgggacctctgtaagacctgggtgcccttcctctctgggcatcctcaccgtc
agcctggtcttaacctccgtga