

IFI6-interferon, alpha-inducible protein 6 Gene

Reference: PTXBC011601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM14C; G1P3; IFI-6-16; IFI616; interferon alpha-inducible protein 6; interferon, alpha-inducible protein clone IFI-6-16; interferon-induced protein 6-16; interferon alpha inducible protein 6		

For research use only.

GENE DESCRIPTION	interferon, alpha-inducible protein 6
GENE ID	2537
ORIGIN SPECIES	Human
ACCESSIONS	BC011601
VECTOR	pDONR223
NCBI SYMBOL	IFI6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggcagaagcggtatcgctttctgtctacgtgctcttcaacttcagtggtgggagggcaggtgaagaaagtgtctggagagctcggacagcggctcgggttctggaaggccctgacctcatggcgtcggaggaggactcgcagtcgccgggctgcccgcgtggg
cttcacggcgccggcatcgccgcaactcgggtggctgctcgtgatgagctggtctgcgatctgaatggggcgccgtgcccgcgggggctagtggccacgctgcagagcctcggggctggtggcagcagcgtcgtcataggtaatattgggtccctgatgggctacgccac
ccacaagtatctgatagtgaggaggtgaggagtag