



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CCND3-cyclin D3 Gene

Reference: PTXBC011616

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

G1/S-specific cyclin-D3; D3-type cyclin

For research use only.

GENE DESCRIPTION	cyclin D3
GENE ID	896
ORIGIN SPECIES	Human
ACCESSIONS	BC011616
VECTOR	pDONR223
NCBI SYMBOL	CCND3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]