

TRIM11-tripartite motif-containing 11 Gene

Reference: PTXBC011629

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase TRIM11; BIA1; RNF92; RING finger protein 92; tripartite motif-containing protein 11; tripartite motif containing 11		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 11
GENE ID	81559
ORIGIN SPECIES	Human
ACCESSIONS	BC011629
VECTOR	pDONR223
NCBI SYMBOL	TRIM11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgaggaccgtgtgcagggtcccgggactggtagagacactgcggaggtttcgaggggacgtgaccttgaccgggacacccgcaacctgagctgacctgtctgaagacaggcggagcgtgcagcggggggacacacggcagccctgccggacagccagagcg
ctttgaccccgcccttgcgtgctggccagagcgttcacctcaggccgacactactggaggtggaggttggggaccgcaccagctgggcccctgggggtgtgcagggaacgtgaacaggagagagagggcgagctgtccgcccgaacgctctggaatcgtgtcttcc
tggggagctattacaattctcgggaacgggccttggctccactccgggacccacccaggcgcgtggggatcttctggaactacgagcgtggacatctcttctacagtgccaccgatgggtcactgctattcatcttcccagatccccttctcggggacgctgcggccctcttcaacc
ctgtccagcagcccgacccgatgactatctgccgccgaaagggtgggtccgggacacccctggctcccagtg