

C22orf29-chromosome 22 open reading frame 29 Gene

Reference: PTXBC011679

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BOP; protein Bop; BH3-only protein; chromosome 22 open reading frame 29		
For research use only.			

GENE DESCRIPTION	chromosome 22 open reading frame 29
GENE ID	79680
ORIGIN SPECIES	Human
ACCESSIONS	BC011679
VECTOR	pDONR223
NCBI SYMBOL	C22orf29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcagaagacacagctagtggcacttggcggttaaacctgcagaaaggggtccctagctgggcacatgccagctcccgacccacagggtggacttctgtgggtaccaggtcagaccagggcaccttggctccctgggtactggaccgttttggccagct
ggcgattacatgtcttccacttggcactatcaggacaacatcagccgtgtctgcgagatctcagggcctaacaggccgagccagggcctggcgagccctacccttggggacctgcccctgctgacgattacgagctcttctgccagatctcaaggagttgttcaagacc
cgaacagttttctgagtagctgtgtgttacctgtccctgcccctggcctccagccagctgccagtgccctcagctgctgtgtgtgaggaatacttagctaggttcttagagggcctggcactcgacatgggtactgccccagggtctttaccagccgcatggccacctgctgt
gtctgggtccaactctgtatctagaagtgtctgttcgagcagcgtgaccaaggagagcaccctgggcccgaaggagccccagctcctgccagttctacatgtagctccaagcctgtgtctgtggaaccagcctcttccagccagaggagcagccccacacctgtccctagact
gtcggagtcagctaactcctgtccagagaccagaccagctcaccagggtccaaaacccagaaaacagaggaggaggttttggagacagaggagaccagagggtgtccttaggtacccacaggaggtgtgtggaggcccccggagacccagagagccaccatttct
cctgggttctga