

C6orf108-chromosome 6 open reading frame 108 Gene

Reference: PTXBC011683

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf108; RCL; dJ330M21.3; 2'-deoxynucleoside 5'-phosphate N-hydrolase 1; c-Myc-responsive protein RCL; deoxyribonucleoside 5'-monophosphate N-glycosidase		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 108
GENE ID	10591
ORIGIN SPECIES	Human
ACCESSIONS	BC011683
VECTOR	pDONR223
NCBI SYMBOL	C6orf108
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtgcatgggtccggggcgcagcgagctgggagcggggagcctggccggccctgtactctgcgggagcattcgcgggacgcgagcagcgtgtacgagcggatcgtgtctcggtgcggcgattcgggacagtgctcaccgagcagctggcg
gcccgcgagctggcgcgcggggaagagctgtgtgggggtgacaggtcatccatgagcaggacctggagtgctgcagcagcgggacgtggtcgtggcagaagtgcacagccatccttgggtgtaggctatgagctggccgggcccgtggccttaacagcggatcctgt
gcctgttccggccgagctgtggcgctgttccgcatgatccggggagcagatggctctcgggtccaggtgtgggactatgaggaggagaggtggaggccctgtggatcgatactcgaggtgatcctcagggcaggtggctgctccctgacccaaccactga