

## CORO2A-coronin, actin binding protein, 2A Gene

Reference: PTXBC011690

|                |                                                                                                                                                                               |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                           | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                          | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CLIPINB; IR10; WDR2; coronin-2A; WD protein IR10; WD repeat-containing protein 2; WD-repeat protein 2; coronin, actin binding protein, 2A; coronin-like protein B; coronin 2A |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | coronin, actin binding protein, 2A                                                               |
| <b>GENE ID</b>          | 7464                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC011690                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | CORO2A                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtcatggcaccacccagtagccgagctccaagttccgtcatgtcttggcaaacccagcagcaagagaaactgctacgactccgtgctatcacccgcagcgttcacgacaaccacttctgtgcccgtgaaccccaactcattgacgtgtgactgagtgctggtggaggggcttctc  
gtcatccctcgaccagacaggggaagttggaccacccactacccaaagtctcgggacagaggcaacgttttggatgtaagtggaaccccttggatgatttggatgcgctcctgttctgaagatgccacaattaagatctggagcatcccaagcagctgtgaccaggaactcac  
ggcctacaggaaggaactctgtggccacgcgcgagtaggctgtgtggagtgccacccacggccgccaacatctcttcagtgctgctatgactacaagtgatgatctggaacctggatacaaggagctgtctatcacaagcccatgagtacaattagctgtcacaagatgt  
gatectctcatgtcttcaacaccaacggcagcctgttggccaccacctgcaagaccgcaagattcgggttattgaccccccagcagggaccgtctcctcaggaggccagctacaaggccacccggccagcaaaagtgtgttctggggaacctgaagaagctgattccacagggc  
acatcccgatggaacaacggcaggtggcctgttggaccaggataacctctgtgctctgatggaggaggacctggacggctcctcgggctgtgtgttccctctatgacgcggacaccagcatgctctacgtgtggggaaggagatggcaacatccgctactacgagtgag  
cgccgacaagcctacctgagctacgtgactgagctaccgtctataaccacagaaggggatcgggtgtatgccaaagagaggactcgacgtgtcctcctgcgagatctccgcttctacaagctgatacaacaaaagcctcatcgagcccatctccatgattgtgccccggcggtca  
gaatcctaccaagaggacatataacctccaacagcagggggccagccctccctgacggccaggagtggtcagcgggatgaatcagagaccaatcctgtgtcccttaggcctggcctgtgagctgtgtgagaccccccactgctgcagagagacatcttcaattccatggcccc  
agcctaccccggtcttgaatcagacagaaaagctgctgcagaagatggctggaggtctctcctcctgttggaggagaagatccaaggtgggcagcagaacacagggctggaggagaagaaaacctggtgacaaatgctgttgacgttttccgaatcccccccaaacagacagag  
aatgagttgctgcagatgttctaccggcaacagaggagatccgaagctccggagctgttgacccagcgagaggtccagggccaacagttggaactggagatcaaaactgcggatgggctcagagcagctctga