



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM43-transmembrane protein 43 Gene

Reference: PTXBC011719

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARVC5; ARVD5; EDMD7; LUMA; transmembrane protein 43		
For research use only.			

GENE DESCRIPTION	transmembrane protein 43
GENE ID	79188
ORIGIN SPECIES	Human
ACCESSIONS	BC011719
VECTOR	pDONR223
NCBI SYMBOL	TMEM43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]