

RAB3A-RAB3A, member RAS oncogene family Gene

Reference: PTXBC011782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB3A, member RAS oncogene family; RAS-associated protein RAB3A; ras-related protein Rab-3A		

For research use only.

GENE DESCRIPTION	RAB3A, member RAS oncogene family
GENE ID	5864
ORIGIN SPECIES	Human
ACCESSIONS	BC011782
VECTOR	pDONR223
NCBI SYMBOL	RAB3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatcgccacagactcgctatggcagaaggagtcctggatcagaacttcgactacatgttcaagattctcatcctggcaacagcagcgtgggcaagacgtccttcttcctgctatgctgacgactgttcacgctgccttcgctcagcaccgtgggcatcgactcaagtc
aagaccatctatcgcaacgacaagagatcaagctgcatctgggacacagcaggcgaagagcgttacccgaccatcaccacgcatactaccggcgctatgggcttcctcatgtatgacatcaccaacgaggaaatcctcaatgcagtgacgagctggccaccagatcaa
gacctactatgggacaatgcccaggtgctgctgtaggaacaagtgtgacatggaggatgagcgggtgtgtcatcagaacgtggccggcagctagctgaccaccttgggttcagttcttggaggcaagcgccaaggacaacattaacgtcaagcagaccttgagcgcctgggtgg
atgtcatctgcgagaagatgtccgagtcgttgacacggcggacccctgcgttcacagcgcccaagcagggcccacagctcagtgaccagcaggtgccaccgcaccagcagctgcgcctgtga