



Reference: PTXBC011819

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent RNA helicase DDX3X; CAP-Rf; DBX; DDX14; HLP2; MRX102; DEAD (Asp-Glu-Ala-Asp) box helicase 3, X-linked; DEAD (Asp-Glu-Ala-Asp) box polypeptide 3, X-linked; DEAD box protein 3, X-chromosomal; DEAD box, X isoform; DEAD/H (Asp-Glu-Ala-Asp/His) b		

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 3, X-linked
GENE ID	1654
ORIGIN SPECIES	Human
ACCESSIONS	BC011819
VECTOR	pDONR223
NCBI SYMBOL	DDX3X
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]