

TSR2-TSR2, 20S rRNA accumulation, homolog (S. cerevisiae) Gene

Reference: PTXBC011825

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSR2, ribosome maturation factor; TSR2, 20S rRNA accumulation, homolog; pre-rRNA-processing protein TSR2 homolog; DBA14; DT1P1A10; WGG1; WGG motif containing 1		

For research use only.

GENE DESCRIPTION	TSR2, 20S rRNA accumulation, homolog (S. cerevisiae)
GENE ID	90121
ORIGIN SPECIES	Human
ACCESSIONS	BC011825
VECTOR	pDONR223
NCBI SYMBOL	TSR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggctgcagaagatgcgcgagctctttccgggctggggtctgcgcggccctggaggcctggccggccttcagatcgtgtggagaatggctcgggggtgtgcacagccaggagaaggccaagtggctgggggtgcagtggaggattctcatgcgcaatgctgac
ttggagctagatgagtggaagacttccttgagagctgtttgaccaacgagtttgatacagttgtggaagacgggagctgccccaggtgagccagcaactgcagaccatgttccaccactccagaggggtgatggggctgctctgaggagatggcctcctgcatcactcagagaaaat
gcaaggctcacagccactgcacttaagacagctagagagactgatgaggatgaagatgatgtggacagttgtggaagagatggaggtcacagctacgaatgatggggctgctacagatggggctgccccagcctgaaccctctgatccagacgctcagactattaaggaaggatata
gtggaagatggctggaccattgtccggagaaaaaatga