

C4orf43-chromosome 4 open reading frame 43 Gene

Reference: PTXBC011842

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0534 protein C4orf43; C4orf43; translation machinery-associated protein 16; translation machinery associated 16 homolog		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 43
GENE ID	55319
ORIGIN SPECIES	Human
ACCESSIONS	BC011842
VECTOR	pDONR223
NCBI SYMBOL	C4orf43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagcaccaagggaagagtcagacggggaagaaagtcacatccatagtagaagcagctcaattacgagagggccacaacaagaaaaaggaattgaagaatgaaaggccttgcgtctcaacctgttggtgaaaaactgcaatggttcaaa
tcaccttgatccccaaaaagagatattcaagaagatgctgtgaactaattgaaggtacttaaatcgattcagcagtgagctggagcagattgagttacataacagtatcaggacagggcaggggagggcgcactgttccgggagaccgtcatcaagcagacgatggagcggga
gagacagcagtttgaggatattggccttgagattccagacattctaaatgcaagtaactggaacatttaggggaatgggactttagctgaagaattaccaacattaaatgagaaaaattgcgctaattgatgcaattccaagacgtgcaagaggaaaactataactgtgaccaag
attgggggaattggaactaacgatgaatcaagtgttcagatgaggaaatgactgcagtgccctaa