



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF343-zinc finger protein 343 Gene

Reference: PTXBC011862

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

dJ734P14.5; zinc finger protein 343

For research use only.

GENE DESCRIPTION	zinc finger protein 343
GENE ID	79175
ORIGIN SPECIES	Human
ACCESSIONS	BC011862
VECTOR	pDONR223
NCBI SYMBOL	ZNF343
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr