



Reference: PTXBC011893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MARVD3; MRVLDC3; MARVEL domain-containing protein 3; MARVEL (membrane-associating) domain containing 3; MARVEL domain containing 3		

For research use only.

GENE DESCRIPTION	MARVEL domain containing 3
GENE ID	91862
ORIGIN SPECIES	Human
ACCESSIONS	BC011893
VECTOR	pDONR223
NCBI SYMBOL	MARVELD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]