

NDUFA4L2-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4-like 2 Gene

Reference: PTXBC011910

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NUOMS; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4-like 2; NADH:ubiquinone oxidoreductase MLRQ subunit homolog; NDUFA4, mitochondrial complex associated like 2		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4-like 2
GENE ID	56901
ORIGIN SPECIES	Human
ACCESSIONS	BC011910
VECTOR	pDONR223
NCBI SYMBOL	NDUFA4L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggagccagcttggggcccgcttaccggcagatcaaaagacatccggggatcatcccgatgatcggcttaactgcctgggcagcgctgcgttacttgcgcgactcgccctcgagccccgacgtctgctgggacagaaagaacaacccggagccctggaacgcctgagcccaatgaccaatacaagttccttgcagtttccactgactataagaagctgaagaaggaccggccagacttctaa