



Reference: PTXBC011913

For research use only.

GENE DESCRIPTION	purinergic receptor P2X, ligand-gated ion channel, 7
GENE ID	5027
ORIGIN SPECIES	Human
ACCESSIONS	BC011913
VECTOR	pDONR223
NCBI SYMBOL	P2RX7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]