

MAGEH1-melanoma antigen family H, 1 Gene

Reference: PTXBC011954

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APR-1; APR1; MAGEH; melanoma-associated antigen H1; apoptosis-related protein 1; melanoma antigen family H, 1; melanoma antigen family H1; restin; MAGE family member H1		

For research use only.

GENE DESCRIPTION	melanoma antigen family H, 1
GENE ID	28986
ORIGIN SPECIES	Human
ACCESSIONS	BC011954
VECTOR	pDONR223
NCBI SYMBOL	MAGEH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctcggggacgaaagagtcggcgccgcgtaatgcgagagccgcagagaagagaaccgcaacaatcgcaaaatccaggcctcagaggcctccgagacccctatggccgctctgtgtagcgagcaccgccgaagacgacctgagcggccccgaggaagaccgagcactcc
agaggaggcctctaccaccctgaagaagcctcgagcactgcccagcacaaaagccttcagtgcccgaggcaattttcagggcaccaaagaaagtcctctgatgtctatattagcgctcatcttcacatcatgggcaacagcgccaaggaagctctggtctggaaagtctggggaagtta
ggaatgcagcctggagctcagcacagcatcttggagatccgaagaagatcgtcacagaagagtttgcgcagaggggtacctgattataaacgggtgccccgtagcagtcgggtggagatgagttcttctggggcccgagcacacgtggaatcgagcaaaactgaaatcatgcatt
ttgtgcgaagggttcgtaaccgatgctctaaagactggccttgaattatgactgggattcggacgatgatgcagaggttgaggctatcctcaattcaggtgctaggggttattccgcccttaa